

UTILITY OF CASPASE-CLEAVED KERATIN 18 TO DIAGNOSE NASH IN PATIENTS WITH OBESITY

Sami Qadri^{1,2}, Hannele Yki-Järvinen^{1,2}, Reda Elkhatab³, Paola Pellegrini³

¹Department of Medicine, University of Helsinki and Helsinki University Hospital, Helsinki, Finland ²Minerva Foundation Institute for Medical Research, Helsinki, Finland ³VLVbio AB, Stockholm, Sweden

Disclosure: SQ & HYJ have no commercial or financial interest in VLVbio

INTRODUCTION

In treating patients with nonalcoholic steatohepatitis (NASH), an ongoing challenge is the paucity of adequately performing noninvasive biomarkers, necessitating the use of a liver biopsy to obtain diagnosis¹. The M30[®] antibody detecting caspase-cleaved keratin 18 (cck18) has emerged as a promising circulating biomarker of liver injury, showing utility to identify patients with NASH². In addition, a cck18-incorporating composite score, MACK-3 (cck18, AST and HoMa), was specifically developed to predict fibrotic NASH^{3,4}, which is the primary inclusion criterion in ongoing clinical trials studying pharmacotherapies for NASH^{3,4}. Although NASH and its complications are especially prevalent in patients with severe obesity and type 2 diabetes⁵, few data exist on biomarker performance within this demographic. Here, we studied the applicability of cck18 circulating fragments to diagnose NASH or fibrotic NASH in a representative obese cohort.

METHODS

We recruited 354 patients aged 25 to 71 years undergoing a liver biopsy during laparoscopic bariatric surgery at the Helsinki University Hospital (Helsinki, Finland). Patients with primary liver diseases other than nonalcoholic fatty liver disease (NAFLD) or significant alcohol consumption (males >30 g/day; females >20 g/day) were excluded. A week before the liver biopsy, the patients underwent clinical examination and blood sampling. Plasma concentrations of cck18 fragments were determined using the M30[®] antibody (M30 Apoptosense[®] cck18 kit [ELISA]; VLVbio, Nacka, Sweden). We calculated MACK-3, the Fibrosis-4 index (FIB-4), and the NAFLD Fibrosis Score (NFS) using their respective formulae^{4,6,7}. Histopathological features of NAFLD were assessed using the NASH Clinical Research Network grading and definitions⁸. Fibrotic NASH was defined as NASH with a NAFLD Activity Score (NAS) \geq 4 and concomitant \geq F2 fibrosis stage³.

RESULTS

Patient characteristics

The mean age of the patients was 51 \pm 10 years and the mean BMI 40.4 \pm 7.2 kg/m². Type 2 diabetes was diagnosed in 43% of the patients. In liver histology, NASH was present in 15%, fibrotic NASH in 6%, and advanced liver fibrosis (stage F3–F4) in 7% of the patients.

Associations between cck18, liver histology, and liver enzymes

Concentrations of cck18 correlated significantly with the full histological spectrum of NAFLD, including steatosis, ballooning, inflammation, fibrosis, and NAS (Figure 1). A linear and positive association was present between cck18 and NAS, but median concentrations decreased at fibrosis stage F4 (Figure 1). In a multiple linear regression analysis, cck18 lost its association with fibrosis stage after adjusting for NAS, confirming its role primarily as a biomarker of liver injury associated with NASH (data not shown). Moreover, cck18 correlated significantly with plasma ALT ($r_s=0.63$, $P<0.001$) and AST ($r_s=0.60$, $P<0.001$).




Figure 1. Associations between concentrations of cck18 M30 and liver histology. * $P<0.05$; ** $P<0.01$; *** $P<0.001$; **** $P<0.0001$.

ccK18 as a predictor of NASH and fibrotic NASH

The ccK18-incorporating MACK-3 score had the highest area under the receiver operating characteristic (AUROC) for NASH (0.84) and fibrotic NASH (0.88) (Table 1, Figure 2). The ccK18-based biomarkers (M30® and MACK-3) had significantly higher AUROCs to diagnose NASH as compared to FIB-4, NFS, ALT, and AST ($P < 0.05$ for all comparisons, DeLong's test), and a significantly higher AUROC to diagnose fibrotic NASH as compared to NFS and ALT ($P < 0.05$).


Table 1. Diagnostic accuracy of the biomarkers (n=354).

Biomarker	NASH	Fibrotic NASH*
	AUROC (95% CI)	AUROC (95% CI)
ccK18	0.84 (0.79–0.89)	0.88 (0.83–0.93)
MACK-3	0.87 (0.82–0.92)	0.89 (0.84–0.95)
FIB-4	0.67 (0.59–0.75)	0.80 (0.70–0.89)
NFS	0.59 (0.50–0.67)	0.63 (0.52–0.74)
ALT	0.73 (0.65–0.80)	0.77 (0.66–0.87)
AST	0.79 (0.72–0.85)	0.83 (0.73–0.92)

AUROC, area under the receiver operating characteristic; CI, confidence interval. *Defined as NASH + NAS \geq 4 + F \geq 2.

CONCLUSIONS

Plasma concentrations of ccK18 associate with the full histological spectrum of NASH. In patients with severe obesity, ccK18 and the related MACK-3 score have an adequate discriminatory ability to identify NASH or fibrotic NASH.

**Figure 2.** Receiver operating characteristic curves for the studied biomarkers to identify (A) NASH or (B) Fibrotic NASH, defined as NASH + NAS \geq 4 + F \geq 2.**REFERENCES**

- Chalasani N, Younossi Z, Lavine JE, et al. The diagnosis and management of nonalcoholic fatty liver disease: Practice guidance from the American Association for the Study of Liver Diseases. *Hepatology* 2018;67(1):328-57. doi: 10.1002/hep.29367 [published Online First: 2017/07/18]
- Lee J, Vali Y, Boursier J, et al. Accuracy of cytokeratin 18 (M30 and M65) in detecting non-alcoholic steatohepatitis and fibrosis: A systematic review and meta-analysis. *PLoS One* 2020;15(9):e0238717-e17. doi: 10.1371/journal.pone.0238717
- Sanyal AJ, Brunt EM, Kleiner DE, et al. Endpoints and clinical trial design for nonalcoholic steatohepatitis. *Hepatology* 2011;54(1):344-53. doi: 10.1002/hep.24376
- Boursier J, Anty R, Vonghia L, et al. Screening for therapeutic trials and treatment indication in clinical practice: MACK-3, a new blood test for the diagnosis of fibrotic NASH. *Aliment Pharmacol Ther* 2018;47(10):1387-96. doi: 10.1111/apt.14621
- Younossi ZM, Golabi P, de Avila L, et al. The global epidemiology of NAFLD and NASH in patients with type 2 diabetes: A systematic review and meta-analysis. *J Hepatol* 2019;71(4):793-801. doi: 10.1016/j.jhep.2019.06.021 [published Online First: 2019/07/08]
- Angulo P, Hui JM, Marchesini G, et al. The NAFLD fibrosis score: A noninvasive system that identifies liver fibrosis in patients with NAFLD. *Hepatology* 2007;45(4):846-54. doi: 10.1002/hep.21496
- Sterling RK, Lissen E, Clumeck N, et al. Development of a simple noninvasive index to predict significant fibrosis in patients with HIV/HCV coinfection. *Hepatology* 2006;43(6):1317-25. doi: 10.1002/hep.21178
- Kleiner DE, Brunt EM, Van Natta M, et al. Design and validation of a histological scoring system for nonalcoholic fatty liver disease. *Hepatology* 2005;41(6):1313-21. doi: 10.1002/hep.20701 [published Online First: 2005/05/26]